

AUTO REFRACTOR

ACCUREF-K

9001

Operations Manual

SHIN-NIPPON

RB-400-B02D

INTRODUCTION

This manual contains information on correct handling and operational procedures as well as safety consideration pertinent to ACCUREF-K 9001.

Before carrying out measurement and/or adjustment, read the instructions thoroughly so that effective operation is ensured. As this constitutes an important reference and user guide, keep it on hand at all times.

NOTE

- The information contained in this manual is subject to change without notice.
- While reasonable efforts have been made in the preparation of this document to ensure its accuracy, you should contact your local distributor immediately, if any quarries arise due to editorial errors or omissions etc.

SAFETY CONSIDERATION

ACCUREF-K 9001 is a Class I, Type B medical instrument as well as LED Class 2 product.






This instrument complies with Medical Device Directive 93/42/EEC.

A great deal of consideration has gone into the design and manufacturing of this instrument with regard to its operational ease, the patient's safety and well-being as well as to the reliability of the product.

For safer and more effective use, however, follow the points described in this manual.

This instrument is designed for professional use.

~ General Definitions of Safety Symbols ~

 CAUTION	General warning. Caution. Risk of danger.		Caution of invisible LED radiation. Avoid exposure to beam. Indicated inside of the device.
	No specific definition given. Denotes general ban or prohibition.		This device is a LED class 2 product.
	General mandatory action.		



CAUTION

- Always take great care when operating ACCUREF-K. Malfunction or damage to the instrument could occur.
- Cut the power immediately if malfunction occurs during operation. Damage to the equipment or personal injury will result. Consult your dealer, if repair work needs to be carried out.



- At no time attempt to remodel or disassemble ACCUREF-K. Damage to the instrument or personal injury will result.
- As ACCUREF-K is a precision optical instrument, operations must be carried out at all times by experienced, authorized personnel. Damage to the equipment or personal injury will result.
- Avoid installation near TV or radio. The reception can be disturbed by electrical noise. Follow the manual for the proper installation.
- Never remove the plug from the outlet if your hands are wet. Electric shock or personal injury could result.
- Make sure the power cord is not damaged. Fire or electric shock may occur.
- Do not touch the optical parts. Measurement accuracy will be adversely affected.



- The power cord must be firmly connected to an electrical ground (safety ground) at the power outlet. Personal injury may result from electric shock, etc.
- If the instrument fails to work properly, you should not try to repair the fault. Consult your dealer immediately.
- The instruction in this manual ensures correct operations.
- Observe the following environmental conditions for used and storage. Avoid dew condensation at all time.

	Temperature	Relative humidity
Use	+10°C to +40°C	30% to 85%
Storage	-10°C to +60°C	below 70%



Avoid the following conditions for storage and use of the instrument.


- Where noxious gases or air pollutants are present.
- Where dust and grit may occur.
- Where oil fumes or greasy substances are emitted.
- Where there are atmospheric concentrations of salt.
- Near gas generation areas and places where dust accumulates.
- Keep in a secure, stable situation. Do not expose to strong vibrations (areas of seismic activity) and sudden shocks (this includes transportation) etc.
- Where there is an inclination of more than 10 degrees.
- Where voltage from the power sources rises or falls sharply during loading.
- Where fluctuations in the voltage of the power source occurs.
- Direct contact with sunlight.

If the instructions above are not followed, damage to the equipment or personal injury will ensue.


CONTENTS

Introduction	1
Safety Consideration	1
Accessories	4
1. Parts Identification	5
2. Conveyance and Handling Procedure	6
3. Installation Environment	6
4. Safeguard Summary	7
5. Preparation	8
5.1 Setting	8
5.2 Applying Power	9
5.3 Standby	10
5.4 Switch Function	11
6. Measurement	12
6.1 Measurement Flow	12
6.2 Alignment	13
6.2.1 R/K, K, and P•K mode	13
6.2.2 R mode	15
6.3 Measurement Results	15
6.4 Print Out	16
6.5 Kerato-Peripheral Measurement	17
6.6 IOL Measurement Function	19
6.7 Menu Screen Setting	20
6.7.1 Each Item Description	21
6.7.2 Optional Functions	24
6.8 Auto Start Function	28
6.9 Data Screen Function	29
6.10 Power Saving Function	30
6.11 Output Terminal	30
7. Tips for Effective Measurement	31
8. Error Messages	32
9. Contact Lens: Base Curve Measurement	33
10. Storage and Maintenance	34
10.1 Reloading Printer Paper	34
10.2 Fuse Replacement	35
10.3 Storage	35
10.4 Confirmation of Measurement Accuracy	36
11. Specification	37


Accessories




Model Eye: 1
 (with a contact lens holder)
 A sticker indicating the
 diopter value is affixed.




Power Cord: 1
 (2.5m)




Printer Paper: 3
 (width 57mm)
 (Two packed and one installed
 into the body)




Fuse: 2
 (T2A 250V)




Pack of Chin Rest Liners: 1
 (1,000 sheets)



Chin Rest Liner Pin: 2



Dustproof Cover: 1



Operations Manual: 1




Use accessories specified by us to avoid any malfunction or failure.



- Extra care should be taken for storage of a model eye. Avoid where the lens of the model eye may be damaged as well as any dusty or humid/steamy environments.
- Avoid direct sunlight, humidity and high temperature when storing printer paper which is a thermal paper.

1. Parts Identification



2. Conveyance and Handling Procedure


When transporting the instrument, make sure that the body has been securely locked.

Center the body onto the base so that their edges are aligned. Push up the lock while turning right until both body and base are firmly fixed to each other.

3. Installation Environment

- 1) Do not expose the instrument's view window directly to the sunlight or bright lighting from other sources.


Great care should be taken and avoid that the examinee is exposed to strong light or glare. The pupil will contract too much for measurement to be carried out.



- 2) Do not operate in places where either dust or rubbish accumulates. Environment with extremes in heat and humidity should also be avoided.

Always follow the environmental requirements below for installation.


Temperature	+10 °C ~ - 40°C
Relative Humidity	30 % ~ 85 %



- 3) Keep away from inflammable or explosive gases as well as storage areas housing medical supplies and chemicals.

- 4) Avoid installing where dew condensation may accumulate. Also, avoid where the radical temperature changes may occur.

- 5) Keep away from sites that may experience strong vibrations or sudden shocks.



- 6) Malfunction is likely to occur if the instrument is improperly stabilized or accidentally overturns. To prevent internal/external damage caused by sudden impact, set the instrument on a solid and secure surface. Do not store in high, 'out of reach' places.

4. Safeguard Summary

1. ACCUREF-K is a precision optical instrument. Always handle with care and avoid dropping it accidentally.
2. Ensure that the instrument is properly grounded when connected to the power source.
3. Do not touch the optical parts with fingers and be sure to avoid dust, as their measuring accuracy could be adversely affected and incorrect values may result.




When dust or fingerprints appear on the optical part, use a soft cloth to wipe off the build-up. In case that the build-up is hard to remove, absolute alcohol is recommended. Take great care when cleaning these parts as they are particularly sensitive and fragile.

4. If the surfaces of the measuring unit and main unit including the control panel are dirty, gently wipe with a dry cloth. For hard to remove stains, a damp cloth or neutral cleanser is recommended.




Avoid using organic solvents that will damage the water based paint finish of the instrument.




5. Clean the chin rest and head rest with the neutral cleanser. For disinfecting them, especially where the examinee may contact, hydrogen peroxide (Oxydol) is recommended.

6. If the instrument is not used for any length of time, remove the power cord from the outlet.

7. When not in use, protect the instrument with a supplied dustproof cover.



8. When the instrument fails to function properly, **never attempt to perform internal service or adjustment.** Contact your nearest registered agent, distributor or retail outlet.




5. Preparation

5.1 Setting

- (1) Set a roll of printer paper in the printer. Refer to '10.1 Reloading Printer Paper' for the procedure.
- (2) Set and fix the chin rest liners with the chin rest liner pins on the chin rest. Refer to the figure on the right.




For sanitary consideration, disposing a sheet of the used chin rest liners after every measurement is recommended.




- Always use the chin rest liners following above.
- For sanitary consideration, disinfecting the chin rest with Oxydol is recommended.

5.2 Applying Power




- (1) Confirm that the power is 'OFF' ().



- (2) Insert the power cord into the instrument's power plug connector. Then insert the plug into a general purpose outlet.




Always make sure that the cable is grounded.



- (3) Turn the power switch 'ON' ().

5.3 Standby

When the power is turned on, the screen as shown below appears on the LCD monitor, which is ready for take measurements.





Icon Description

(Normal Measurement)

Icon	Description
	Increment of refraction data
	Indication of corneal vertex distance Options are 0, 10, 12, 13.5 and 15mm.
	Mark of cylindrical value
	Indicates an eye to be measured.
	Refractive data display area. Values of S, C and A are indicated.
	Keratometry data display area. Indicates values of R1, R2, A. It is also available to indicate the results as K1, K2, and KC, setting in the mode selection screen.
	IOL measurement mode indication
	Selects measurement mode: refraction & keratometry continuous measurement (R/K mode), refractive measurement (R mode), keratometry measurement (K mode), and kerato-peripheral measurement (P.K. mode).
	Switches over to Menu screen.
	Displays the measurement result on the screen and prints it out.

5.4 Switch Function

Operating switches under a LCD monitor corresponds to the icons displayed on the bottom of the screen. For normal measurement, each switch functions as below.



Operation of Joystick

- Forward: moves the measuring unit toward an examinee.
- Backward: moves the measuring unit toward an examiner.
- Right/Left: moves the measuring unit to the right or left respectively.
- Moves the measuring unit up when rotated to the right and down when to the left.

6. Measurement

6.1 Measurement Flow

- (1) Rotate an anti-sliding lock clockwise to release it.
- (2) Have the examinee place his/her chin on the chin rest and his/her forehead against the head rest.
Then, have him/her see a target.



Uncomfortable posture may cause fatigue to the examinee during measurement.
Vertically adjust the optical table or the chair to avoid that.

- (3) Check from the side and adjust the chin rest until the examinee's eye level is in line with the eye mark.
- (4) When the eye appears on LCD monitor, carry out alignment for correct measurement.



Refer to '6.2 Alignment' for detail procedure.

- (5) Press the measurement start switch when the alignment mark reaches a center of the reticle mark, and take measurement.



ACCUREF-K is equipped with high-speed measurement function.
This function allows you to measure continuously if you keep pressing the measurement start switch during measurement.

6.2 Alignment

For more precise measurement, ACCUREF-K indicates an alignment mark to assist alignment and focus indicators to indicate the focus status on the monitor.

6.2.1 R/K, K, and P•K mode


With this mode, you need to take two steps to measure with more precise work distance: focus detection and focusing. For focus detection, focus indicators are **red**, and then change to **blue** at focusing.

【Focus Detection】


- (1) Match the center of the eye to a reticle mark. When you try to bring the eye into focus, alignment mark (+) will appear.




If the eyelid is over the kerato rings, urge the examinee to open the eye bigger.





- (2) Operate the joystick to bring the alignment mark (+) into the center of the reticle mark.



- (3) Keeping the alignment mark (+) in the center of the reticle mark, move the joystick toward (or away from) you so that the eye image will be defocused. Consider that the focus point of the eye as a starting point at this time. Along with the joystick's movement, the focus indicators (▬▬) will move outward, and then, stop.



- (4) Move the joystick away from you (or toward you) so that the eye will be focused, and then, defocused. At this time, the focus indicators () will move inward, and then, outward. Stop moving the joystick when the indicators move outward and stop. Focus detection will complete when color of the focus indicators change from red to blue.





The joystick can be moved only in the range where the alignment mark (+) remain on the screen.



The color of focus indicators may not change to blue after the procedure above if the examinee blinks too frequently and/or you move the joystick too fast. In such a case, move the joystick slowly back and forth a few times.

【Focusing】

- (1) When the focus indicators () change to blue, move the joystick until the eye image becomes clear. The focus indicators will move inward and disappear once the appropriate position is achieved.
- (2) When the indicators disappear, press the measurement switch to start measurement.
(For auto measurement, disappearance of the indicators will trigger measurement.)



【Indicators to be disappeared】



If you measure with the focus indicators on, the result will not be reliable. Always measure after you confirm that the focus indicators have disappeared. If the indicators reappear after one measurement, you need to realign.





It is essential to achieve correct alignment for more precise measurement. It is strongly recommended to practice and be accustomed to alignment and focusing before measurement.

6.2.2 R mode


Focus indicators will not be indicated on the screen for R mode.

Bring alignment mark (+) into the center of reticle mark.
Focus on the eye and press the measurement switch.

(For auto measurement, measurement will start once the focus on the eye is achieved.)



6.3 Measurement Results



PD result is indicated after both right and left eyes are measured. The order of the eye to be measured is not important.

NPD is indicated when some value is set in 'W-D (cm)' of menu screen only.

6.4 Print Out

Normally you can print out the measurement result after the measurement, using PRINT switch.

For refractive measurement, a maximum of ten data for each eye can be saved and the most reliable value among them is indicated as an optimum value. However, it is printed out only when each eye is measured more than three times.

For keratometry measurement, an optimum value only is printed out.

Format of the printout, ALL, ECONO and OFF, can be set on the menu screen.

- * ALL : in case of refractive measurement, prints out a maximum of ten data for each eye.
For other measurements, respective optimum values only are printed out.
- * ECONO: prints out only optimum values of all measurement data.
- * OFF : prints out no data.

Sample Print Out when PRINT FORM is set to ALL		Sample Print Out when PRINT FORM is set to ECONO	
Message Area	<div> <div>No. 00001</div> <div>NAME 1999 11 30 14:30</div> <div>VD=12</div> <div><R> SPH CYL AX</div> <div>- 4.75 -0.25 62</div> <div>- 4.75 0.00</div> <div>- 4.75 0.00</div> <div>- 4.75 0.00</div> <div><R> mm D AX</div> <div>R1 7.59 44.50 120</div> <div>R2 7.57 44.50 30</div> <div>AVE 7.58 44.50</div> <div>CYL 0.00</div> <div><L> SPH CYL AX</div> <div>- 4.50 -0.75 90</div> <div>- 4.50 -0.75 90</div> <div>- 4.50 -0.75 89</div> <div>- 4.50 -0.75 90</div> <div><L> mm D AX</div> <div>R1 7.59 44.50 120</div> <div>R2 7.57 44.50 30</div> <div>AVE 7.58 44.50</div> <div>CYL 0.00</div> <div>PD = 64</div> <div>SHIN-NIPPON</div> <div>ACCUREF-K 9001</div> </div>	Date and Time Display	
Refractive Data		Right Eye Data	<div> <div>No. 00001</div> <div>NAME 1999 11 30 14:30</div> <div>VD=12</div> <div><R> SPH CYL AX</div> <div>- 4.75 0.00</div> <div><R> mm D AX</div> <div>R1 7.59 44.50 120</div> <div>R2 7.57 44.50 30</div> <div>AVE 7.58 44.50</div> <div>CYL 0.00</div> <div><L> SPH CYL AX</div> <div>- 4.50 -0.75 90</div> <div><L> mm D AX</div> <div>R1 7.59 44.50 120</div> <div>R2 7.57 44.50 30</div> <div>AVE 7.58 44.50</div> <div>CYL 0.00</div> <div>PD = 64</div> <div>SHIN-NIPPON</div> <div>ACCUREF-K 9001</div> </div>
Keratometry Data		Optimum Value (Indicated when each eye is measured more than three times)	
Refractive Data		Left Eye Data	
Keratometry Data			
Pupil Distance			

Message Area


You can print out registered characters in the range of 22 character/line × 2 lines in the message area. For registering characters, refer to '6.7.2 Option Function: 1. Message Input Function'.

6.5 Kerato-Peripheral Measurement

ACCUREF-K has a function to measure not only a center but also peripheries of a cornea.

Operation

- Press a measurement mode selection switch to switch over to P.K. mode.
P.K. measurement screen will appear, and 'measurement guiding mark' (see NOTE below) which indicates the measurement position will appear on the top of the screen.



Measurement position :

Normally measurement is taken in the following order:

C S T I N.

C (Central)	: center
S (Superior)	: top
T (Temporal)	: temple side
I (Inferior)	: bottom
N (Nasal)	: nasal side



Measurement Guiding Mark

A measurement guiding mark changes its color and state to signal the measurement position and the end of measurement. Each color and state means as follows.

No color:	where have not been measured yet.
Blue / Blinking:	where is going to be measured.
Yellow / On:	where measurement is complete.
Blue & Yellow /Blinking:	where is measured once but going to measure again.


- Measure the center, following kerato measurement procedure.

When the measurement is complete, the center of the measurement guiding mark (blue/blinking) will change to yellow illumination. Then, the next measurement position (top) starts blinking.



Always start to measure from the center. However, it is not necessary when the center is already measured in R/K mode or Kerato mode.


Position to be measured (blue/blinking)



- (3) Start to measure the corneal periphery.
Have the examinee look at the fixation target that is turned on. A measurement position of the measurement guiding mark turns blue and blinking at this time. Carry out alignment and press the measurement switch. When the measurement is completed, the blue blinking indication will change to yellow illumination. Then, the next measurement position will start blinking.



For peripheral measurement, bring the kerato rings to the center of the reticle mark, not to the center of pupil.



- (4) Complete the peripheral measurement for all four positions. When it is completed, all positions of the measurement guiding mark will turn yellow and illuminate. Proceed to measure the other eye, following the procedure above.



- When you wish to measure some position again, use P.K. target selection switch to move the cursor to the position you wish to measure.
- When you can not take any data, or you do not need all data of the peripheral measurement, you can skip some position using P.K. target selection switch. In such a case, the result of the position measured only is displayed. The center, however, must be measured all the time.

Printout Sample of P.K. Measurement

Peripheral Measurement Result				<R>	mm	D	AX
				R1	7.59	44.50	120
				R2	7.57	44.50	107
				AVE	7.58	44.50	
				CYL		0.00	
				S	7.84		
				V	7.51		
				T 7.66 - H	7.55	- N 7.71	
				I	7.84		
				V	0.554	H	0.396
				e		AVE	0.475
				<L>	mm	D	
				R1	7.50	44.50	

Data of Corneal Center only
When R/K or Kerato measurement is also taken, only optimum values are printed.

Average of Vertical and Horizontal Eccentricities

Horizontal Eccentricity

Vertical Eccentricity


6.6 IOL Measurement Function

ACCUREF-K has the function to measure the IOL (intraocular lens) implanted eye.

When measuring the IOL implanted eye, press IOL switch on the front panel of the body.

Then, icon for IOL mode is indicated on the left side of the LCD monitor.

IOL measurement mode indication



【Measurement Screen when IOL mode is set】

No. 00001			
NAME			
1999	11	30	14:30
VD=12			
<R>	SPH	CYL	AX
IOL	- 4.75	0.00	
IOL	- 4.50	-0.25	60
IOL	- 4.50	-0.25	131
	- 4.50	-0.25	60

【Print Out Sample】

Measurement results of IOL measurement mode have 'IOL' marks on the left side of each data.

IOL measurement mode will be cancelled:

- when IOL switch is pressed again
- when the measurement mode is switched over to either K mode or P.K. mode
- when print switch is pressed
- when the eye to be measured is switched from right to left or vice versa.
- when the power is turned off.



If IOL mode is not set when IOL implanted eye is measured, error message may appear and measurement may fail.

6.7 Menu Screen Setting

Standard measurement mode is preset to be ready to use. However, you can easily alter the setting if you wish.

To enter the menu screen, press a menu switch below the LCD monitor.

【Menu Screen】

STEP	0.25	0.12
VD(mm)	0 10 12 13.5 15	
CYL	- + ±	
START	MANUAL	AUTO
REF	NORMAL	QUICK(3)
KERATO	RADIUS	DIOPT
PRINT FORM	ALL	ECONO OFF
PRINT OUT	TOP	BOTTOM
DATA SCREEN	ON	OFF
W-D(cm)	OFF 30 35 40 45	
TARGET LIGHT	BRIGHT	NORMAL DARK
SCREEN ADJ.		

the First Screen





SAVE(min.)	OFF 3 5 10
BUZZER	HIGH LOW OFF
OPTION	MESSAGE No. RS232C
DATE FORM	YMD DMY MDY
DATE	1999/11/30
TIME	14:30:22

the Second Screen

Change of Switch Function


Each switch will change its function when you enter the menu screen.

Follow the icons indicated on the bottom of the screen, which corresponds to each switch (see below).

- ◆  Moves the cursor downward on each setting menu.
- ◆  Moves the cursor upward on each setting menu.
- ◆  Selects the item in each setting menu. The cursor moves to the right.
- ◆  Completes the setting and returns to the measurement screen.

6.7.1 Each Item Description

【The First Screen】

- **STEP** Selects the increment for refractive measurement.
 - **VD(mm)** Selects corneal vertex distance.
 - **CYL** Selects the sign for astigmatism.
 - **START** Selects how to start measurement. Refer to '6.8 Auto Start Function' for detail.
 - MANUAL**: begins measurement every time a measurement switch is pressed.
 - AUTO** : begins measurement automatically when measurement requirement is met. Three measurements for each eye are taken, and the data is automatically printed out with completion.
 - **REF** Selects refractive measurement method.
 - NORMAL**: takes one measurement when the measurement start switch is pressed one time.
 - QUICK** : takes measurements continuously for the number of times preset when the measurement start switch is pressed one time. A maximum of ten measurements can be taken.
- 

NOTE REF setting is invalid when START is set to AUTO.
- **KERATO** Selects the unit for radius corneal curvature.
 - RADIUS** : Radius of corneal curvature
 - DIOPT** : Corneal refraction
 - **PRINT FORM** Selects printout format. Refer to '6.4 Print Out' for detail.
 - ALL** : prints all data out.
(maximum ten data each in case of refractive measurement.)
 - ECONO** : prints only optimal values.
 - OFF** : no print outs.
 - **PRINT OUT** Selects direction of data printed out.
 - **DATA SCREEN** Displays measurement results saved in the memory on the screen. Refer to '6.9 Data Screen Function' for detail.
 - ON** : displays all measurement results on the screen.
 - OFF** : do not display the result on the screen

- **W-D (cm)** Sets near work distance. When measurement is taken with this item set, near pupil distance is automatically computed to indicate it on the screen and the printout.





This function becomes valid when refractive measurement is taken.



- **TARGET LIGHT** Sets brightness of the target.
 BRIGHT : brighten the target.
 NORMAL: normal setting.
 DARK : darken the target.

- **SCREEN ADJ.** Adjust or alter brightness of LCD monitor.



Switch Functions will change in this item as below.

-  (R·K>R>K> P.K. Switch) : makes the monitor brighter
-  (IOL Switch) : makes the monitor darker



When you finish adjusting or altering, move the cursor to any other items except for "SCREEN ADJ.", using  or  switch.




【The Second Screen】

- **SAVE (min.)** Selects switchover time (in minute) to activate power saving function.
- **BUZZER** Sets volume of buzzer at measurement.
 HIGH : turns volume up.
 LOW : turns volume down.
 OFF : no buzzer
- **OPTION** Selects and sets optional functions. Refer to '6.7.2 Optional Function' for detail.
 MESSAGE : shifts to the screen for registering message.
 No. : shifts to the screen for setting examinee's number.
 RS232C : shifts to the screen to set RS232C transmission parameters.
- **DATE FORM** Selects display form of date.
 YMD: year / month / date
 DMY: data / month / year
 MDY: month / date / year

- Sets and corrects time.

If you wish to correct date or time, move the cursor to the number you wish to change.




When you complete all settings or changes, press  or  switch to move the cursor to any items other than TIME. Then, press  switch to return the measurement mode.

6.7.2 Optional Functions

When you select the function you wish to set from 'OPTION' of menu screen, you can enter each option screen.


【Each Option Screen and Description】

I. Message Input Function



With this function, you can print out registered messages in the area of 22 characters / line × 2 lines.

When you select this menu, the screen on the left will appear, and function of each switch will change as shown on the left.




- (1) With or switch, select the character you wish to input and confirm with switch.

Then, the input cursor moves to the next input area to be ready for the next input.



When you need to change the character already inputted, press switch to move the input cursor to the character you wish to input. Then, you can write over.




〔 Print Out Sample 〕

- (2) When the setting is completed, move the input cursor to 'END' and press switch to go back to the menu screen.

II. No. Function


You can set or change the examinee's number, and select whether the number is displayed on the screen and whether the number is printed out.

When you select this menu, the screen below will appear and function of each switch will change as following.



Resetting Examinee's Number

When you wish to reset the examinee's number, move the cursor to RESET of No.SET and press the measurement start switch.




【Screen when DISPLAY is set to ON】


(1) Move the cursor to the item you wish to set/change with switch and execute with switch.

(2) When you complete the setting, press switch to go back to the menu screen.

【Print out when PRINT is set to ON】



Examinee's Number



III. RS232C Setting Function

BAUDRATE
 CHARACTER
 PARITY
 STOP BIT

9600
 8
 NONE
 1

↓

↔

↶

When you select this menu, the screen on the left will appear and function of each switch will change as shown.

Moves the cursor downward to select the item you wish to change

Return switch to go back to the measurement mode.

Changes the setting of the item selected.
 The display changes every time the switch is pressed.
 Press the return switch at any given setting.

With this function, you can send the measurement data to an external computer through an interface. The data is sent using ASCII CODE.

【Each Item Description】

1. BAUDRATE : Selects a kind of transmission rate for a serial interface.

BAUDRATE to be selected	Standard Setting
38400bps	
19200bps	
9600bps	
4800bps	
2400bps	

2. CHARACTER : Selects a kind of data bit for a serial interface

CHARACTER to be selected	Standard Setting
8	
7	

3. PARITY : Selects a kind of data check transmitted from a serial interface.

PARITY to be selected	Standard Setting
EVEN	
ODD	
NONE	

4. STOP BIT : Selects a kind of code to terminate data output from a serial interface.




STOP BIT to be selected	Standard Setting
2	
1	

BAUDRATE	9600
CHARACTER	8
PARITY	NONE
STOP BIT	1

↓


↔

↶

- (1) Move the cursor to the item you wish to set/change with  switch and execute with  switch.
- (3) When the setting is completed, press  switch to go back to the menu screen.

【Connector Pin Layout and Connection】

(1) Connector Pin Layout



(2) Connections

ACCUREF-K		Computer (PC-AT compatible)	
TXD	2	2	RXD
RXD	3	3	TXD
DSR	6	4	DTR
SG (GND)	7	5	GND
DTR	20	6	DSR



Use a shield type of cable for the connection cable in order to protect the output data from noise.

If you have any inquiries about changes in transmission parameters and/or detail connections, please direct them to the relevant agent or distributor.

6.8 Auto Start Function

ACCUREF-K has Auto Start function.

This function starts measuring automatically when alignment meets the measurement requirement, and also prints out automatically when the measurement of both eyes completes.



When Auto Start Function is activated, measurement is always taken continuously for three times regardless of REF setting in the menu screen.

STEP	0.25	0.12
VD(mm)	0 10 12 13.5 15	
CYL	- + ±	
START	MANUAL AUTO	
REF	NORMAL QUICK(3)	
KERATO	RADIUS DIOPT	
PRINT FORM	ALL ECONO OFF	
PRINT OUT	TOP BOTTOM	
DATA SCREEN	ON OFF	
W-D(cm)	OFF 30 35 40 45	
TARGET LIGHT	BRIGHT NORMAL DARK	
SCREEN ADJ.		

- (1) When you wish to measure with auto start function on, set 'START' in the menu screen to 'AUTO.'

- (2) Then, press switch to go back to the measurement screen.

- (3) Auto Start mode is activated when 'AUTO' and arrows below 'AUTO' appear on the top of the screen (see on the left).

You can measure either left or right eye first.

- (4) Carry out alignment. Measurement will automatically start when alignment is achieved.



For alignment, refer to '6.2 Alignment'.

- (5) When the measurement of one eye finishes, an arrow on either left or right side you have measured will disappear.



Arrows and signal the operator completion of measurement by disappearing.

- (6) Take a measurement for another eye following the steps above.

- (7) When measurement of both eyes are over, the arrow will disappear and the result will be automatically printed out.



The data printed out depends on the setting of print out format.

When the data screen function is set on, representative data is displayed on the screen after both eyes are measured.

- (8) If you still wish to measure with auto start function, realignment is necessary.







When you wish to cancel auto start function, change the setting of 'START' back to 'MANUAL' in the menu screen.

6.9 Data Screen Function


Data screen function allows you to check measurement results saved in the memory on the screen.

When displaying measurement data on the screen


STEP	0.25	0.12
VD(mm)	0 10 12	13.5 15
CYL	-	+
START	MANUAL	AUTO
REF	NORMAL	QUICK(3)
KERATO	RADIUS	DIOPT
PRINT FORM	ALL	ECONO OFF
PRINT OUT	TOP	BOTTOM
DATA SCREEN	ON	OFF
W-D(cm)	OFF	30 35 40 45
TARGET LIGHT	BRIGHT	NORMAL DARK
SCREEN ADJ.		

** RIGHT**					
R) SPH	CYL	AX		mm	D
-4.75	-0.25	62	R1)	7.59	44.50
-4.75	0.00		R2)	7.57	44.50
-4.75	0.00		AVE	7.58	44.50
-4.75	0.00		CYL		0.00
-4.75	-0.25	62	S	7.58	T 7.66
-4.75	0.00		V	7.51	H 7.55
-4.75	0.00		I	7.84	N 7.71
-4.75	0.00				
-4.75	0.00		e(ver)	0.554	
-4.75	0.00		e(hol)	0.396	
-4.75	0.00		e(ave)	0.475	



PRESS START SW ->
 RETURN




- (1) Select 'ON' in 'DATA SCREEN' field of the menu screen.



NOTE

When 'DATA SCREEN' is set ON, the setting of 'PRINT FORM' (see 6.7.1 Each Item Description) is invalid.

- (2) After the measurement, press  switch and the data screen on the left will appear.





NOTE

You will see 'I' mark at the left side of each refractive data as shown below when you measure with IOL measurement mode.

Example: SPH CYL AX
 I - 4.75 - 0.25 62

- (3) Displaying the data on the screen, press  switch to change over the display as below.




- (4) When you wish to print out the data, press  switch.

- (5) Press the measurement start switch to go back to the measurement screen.

6.10 Power Saving Function

Power saving function will start operating when switch operation is suspended with the power on.
(The switchover time (in minute) can be selected at 'SAVE' on the mode selection screen.)

To return to the measurement mode, press any switch (any switch on the front panel or the measurement start switch).



6.11 Output Terminal

Video Terminal

This terminal outputs an NTSC video signal.

If you connect an external monitor such as a portable TV or similar unit to ACCUREF-K with pin plug (video) cable, you can observe and check the same image that appears on the internal monitor screen of ACCUREF-K simultaneously.

7. Tips for Effective Measurement

- (1) Do not allow external light to directly penetrate the room.
- (2) Fluctuation of values during measurement may occur if the examinee looks something other than the target. Urge the examinee to concentrate on the target set in front.
- (3) Talk to the examinee in a relaxed and friendly manner, so as to allay any fear or doubt they may have.
- (4) Inappropriate height of a chin rest or a chair will cause the examinee fatigue. Adjust the (optional) instrumental table to establish the most comfortable and convenient position for the examinee.
- (5) When the eyelash or eyelid interfere measurement, error will occur in measurement. Urge the examinee to keep his/her eye wide open.
- (6) Tear residue or eye mucus, etc. trapped on the corneal surface may cause measurement errors. Check the surface with LCD monitor, and if you see something moving when the examinee blinks, remove it before measurement.
- (7) When the pupil of the eye to be examined is smaller than the minimum pupil diameter measurable, correct measurement will be impossible. When the pupil is too small to take correct measurement, make the surroundings (room) or the target darker to allow the pupil to dilate as much as possible.
- (8) If the examinee's head moves during measurement, AXIS value will be adversely affected. Ask him/her to maintain correct posture.

8. Error Messages

ACCUREF-K automatically evaluates measurement condition or result and indicates error messages if it is invalid. Error messages also appear when abnormality is detected in its operational system.

When any error message appears, always check the system with a supplied model eye. If it appears when no abnormality in system is detected, check the measured eye for eye disease or problem.


Message	Cause/Explanation	Corrective Action
RETRY	Failed to capture eye image because the examinee blinks or moves during measurement or the examined eye has an eye disease.	Realign precisely and take measurement again. Consult your dealer immediately if the message reappear. Do not try to repair by yourself.
SPH OVER	Exceeded spherical measurement range (−25D to +25D). (when VD = 0)	
CYL OVER	Exceeded cylindrical measurement range (0 to ± 10 D) .	
Motor fault	Detected abnormality in motor control system.	Cut the power and turn it back on. Consult your dealer immediately if the message reappear. Do not try to repair by yourself.
EEPROM fault	Failed to initialize.	
Print Head Up	Printer head is up.	Close the printer head.
Print Head Heat Over	Printer head is overheating.	Cut the power and stop using until the head cools off. Consult your dealer immediately if the message reappear. Do not try to repair by yourself.
Print Cutter fault	A paper jam occurred at a printer cutter or the printer cutter did not move for some reason.	Always cut the power , and check on the paper jam. If the message reappears without the paper jam, consult your dealer immediately. Do not try to repair by yourself.
Paper Empty	No printer paper.	Set the printer paper. See '10.1 Reloading Printer Paper.'

9. Contact Lens: Base Curve Measurement

You can measure a base curve of a hard contact lens with ACCUREF-K.

To do so, attach a contact lens onto a contact lens holder of the model eye as following.

- (1) Put a small amount of water on a concave side of the contact lens holder.
- (2) Place the contact lens so that its convex side faces the holder.




- (3) Confirm the contact lens is firmly adhered to the holder and does not slip down, set the model eye unit to measure.

10. Storage and Maintenance

10.1 Reloading Printer Paper

- 1) Remove a printer cover and take a printer paper shaft out.




- 2) Lift a printer cutter part and pull up a head-up lever.

- 3) Set a roll of printer paper, paying attention to direction of the paper rolled up.



NOTE

The loose end of the paper should be drawn forward in a counterclockwise direction.




- 4) Insert the loose end of the paper behind a rubber roller. Press print switch more than one second to feed the paper out.




NOTE

Never pull the paper out. Pulling the paper may cause the paper output incorrectly or jammed.




- 5) Draw the paper through the printer cutter slot and lower the head-up lever and the printer cutter back to the original positions.



- 6) Draw the paper a cutter slot of a printer cover and reattach the printer cover to complete the procedure.





Always use the specified printer paper only. The paper other than the specified one may cause a paper jam or fade in printing.



10.2 Fuse Replacement

When a fuse is blown, remove a fuse holder at the side of a main unit for replacement.
Pushing the fuse holder, rotate it in the direction of the arrow below and you can remove it.



Rotate the holder counterclockwise.

10.3 Storage

(1) Points to be checked for long-term storage

- Turn the power switch OFF.
- Remove the power cord from the outlet.
- Lower the optical unit to the bottom (original position).
- Secure the body with the anti-sliding screw lock.
- Put a dustproof cover on the optical unit.

(2) Notes on storage environment

Avoid storage under the following conditions.

- Where dust accumulates.
- Where water may get on the unit.
- Where temperature and humidity are high.
- Where sunlight directly contacts.
- Unstable and/or high place.

Always follow the environment conditions below for storage.

Environmental Conditions	
Temperature	Relative Humidity
- 10 ~ + 60	below 70%



Always check above whenever you store or do not use ACCUREF-K for a long time.
When you reuse the instrument after long-term storage, operate according to instruction in '5. Preparation.'

10.4 Confirmation of Measurement Accuracy


It is extremely important to check operation and accuracy of the instrument using a supplied model eye.

We recommend you to check accuracy periodically.

When the measurement result of the model eye falls anywhere within the tolerance listed below, measurement should be considered reliable and accurate. When the result exceeds the tolerance, contact your dealer immediately.

Model Eye Data		
SPH	CYL	R
Indicated value ± 0.25	0 ± 0.25	Indicated value ± 0.03

Precise value of the supplied model eye is indicated on the model eye stand (VD =12).



NOTE

Note for setting of model eye

- Remove a contact lens holder cap and set a model eye, checking it is not inclined in any directions.

When the model eye is inclined, CYL value can not be correctly measured.

- Set the model eye at the position where an alignment mark is located inside a reticle mark and focus guiding marks disappear from the screen (in other words, good focusing is achieved).
- When conditions above have been satisfied, proceed to measurement.

11. Specification

Refractive Measurement Range	Sphere (S)	-25 ~ +25D	(step: 0.12/0.25D)
	Cylinder (C)	0 ~ ± 10D	(step: 0.12/0.25D)
	Axis (A)	0 ~ 180°	(step: 1°)
Measurement Accuracy	Sphere	between - 10 ~ + 10D: ± 0.25D beyond ± 10 : ± 0.5D	
	Cylinder	± 0.25D	
Keratometry Measurement Range	Radius of Corneal Curvature	5.0 ~ 10.0mm	(step: 0.01mm)
	Corneal Refraction	33.75 ~ 67.5D	(step: 0.12/0.25D) (a refractive index of cornea: n = 1.3375)
	Cylinder	0 ~ ± 9D	
	Axis	0 ~ 180°	(step: 1°)
Vertex Distance	0, 10, 12, 13.5, 15mm		
Minimum Pupil Diameter	2.3 mm		
PD Measurement	Measurement range	85mm	(step: 1mm)
Measurement Time	Refractive measurement	approx. 0.07 sec.	
	Keratometry measurement	approx. 0.07 sec.	
Printer	Thermal line printer with automatic cutter (paper width 57mm)		
Internal Monitor	5.6 inch LCD monitor (color)		
Shifting Range for Sliding Body	back/forth ± 17mm	right/left ± 43mm	up/down 17 mm
Vertical Adj. Range for Chin Rest	± 30mm		
Dimensions	240 mm (W) × 418 mm (D) × 438 mm (H)		
Weight	approx. 15kg		
Data Output	RS232C interface		
	Video Terminal		
Power Source	100 ~ 240V		
	50 / 60Hz		
Consumption	80VA		
Power Saving Function	OFF , 3 , 5 , 10min. (switchable)		

Production Year

The second digit of serial number represents a production year of each instrument. The serial number is indicated on the nameplate that is on the side of body. See below.

5 0 P 0 1 0 3

This number is the last digit of the production year.
In this example, the production year is 2000.

SHIN-NIPPON COMMERCE, INC.

Konwa building
12-22 Tsukiji 1-chome, Chuo-ku
Tokyo 104-0045, JAPAN
TEL:(03)3542-4714 FAX:(03)3545-0507



Distributed by